

FISCHER 426Q THERMO-HYGROGRAPH

Ask About
Traceable Calibration

Temperature

Calibrator

Pressure

Humidity

Recorder &
Data loggers

Air Flow

Level & Flow

Portable
Instruments

EX Proof / IS

Environment

Aerospace



Features

- The carrying-handle for the ease of portability
- Easily detachable housing simply slides out
- Simple removal of the drum
- Attaching recording graph paper is simple and easy
- Sensor of the thermograph is a U-shaped bimetal of high specific deflection

SPECIFICATIONS

Temperature:	-5 to 55°C
Accuracy:	±0.5°C
Humidity:	0 to 100% RH
Accuracy:	±3% RH (20 to 100 % RH)
Rotation (Days):	1, 7, 31 days
Dimension:	290 (W) x 145 (D) x 260 (H) mm
Weight:	3.1 kg
Standard Accessories	1 pack of 7-day chart paper (60 sheet) 2 nos of pens (black) 1 x AA battery Certificate of Conformance (COC)

MODEL

DESCRIPTION

TH426Q	Thermo-Hygrograph
<u>Optional</u> CAL	Traceable Calibration Certificate
R-ME-PEN-B	Cartridge pen for Thermohygrograph (Black)
R-ME-PEN-R	Cartridge pen for Thermohygrograph (Red)